



3233 15th Street South
 Fargo, ND 58104, USA
 fax: 1 (701) 499-7298
 phone: 1 (701) 297-9256
 toll-free: 1 877 PURE-DNA
 email: DNA@aldevron.com

Thank you for your interest in Aldevron. To place an order, fax or mail this form with the necessary information. Orders cannot be started until the appropriate form is completed and the required information is supplied. Incorrect or incomplete information may incur additional charges and prolong turnaround time. Specific order forms are available for download in PDF format from www.aldevron.com/forms, or can be faxed or emailed upon request.

Client / Company: _____

High-Expression Reporter Plasmids

Genlantis has engineered a new series of reporter plasmids based on the gWiz™ high-expression vector system. These plasmids are designed to produce the highest levels of gene expression in a wide range of mammalian cells and tissues. They contain a proprietary modified promoter, followed by intron A from the CMV immediate early gene and a high-efficiency artificial transcription terminator. The expression vectors are constructed in a modified plasmid backbone to achieve enhanced levels of transgene expression in mammalian cells.

Aldevron offers ultra-pure, milligram amounts of these plasmids for use by Gene Therapy and DNA Vaccine researchers. The ready-to-use plasmids are endotoxin-free and can be applied directly to most *in vitro* and *in vivo* applications. Turnaround time is approximately 24 to 48 hours. These plasmids are shipped in water at approximately 5mg/ml.

Catalog #	Description	Size	Price	X Qty	Subtotal
5000	gWiz™ High-Expression Luciferase vector (gWiz™ Luciferase)	5 mg	\$600		
5001	gWiz™ High-Expression Luciferase vector (gWiz™ Luciferase)	10 mg	\$850		
5002	gWiz™ High-Expression beta-galactosidase vector (gWiz™ Beta-gal)	5 mg	\$600		
5003	gWiz™ High-Expression beta-galactosidase vector (gWiz™ Beta-gal)	10 mg	\$850		
5004	gWiz™ High-Expression Secreted Alkaline Phosphatase vector (gWiz™ SEAP)	5 mg	\$600		
5005	gWiz™ High-Expression Secreted Alkaline Phosphatase vector (gWiz™ SEAP)	10 mg	\$850		
5006	gWiz™ High-Expression Green Fluorescent Protein vector (gWiz™ GFP)	5 mg	\$600		
5007	gWiz™ High-Expression Green Fluorescent Protein vector (gWiz™ GFP)	10 mg	\$850		
5008	gWiz™ High-Expression control vector (gWiz™ Blank)	5 mg	\$600		
5009	gWiz™ High-Expression control vector (gWiz™ Blank)	10 mg	\$850		
5010	gWiz™ High-Expression HBsAg vector (gWiz™ HBsAg)	5 mg	\$600		
5011	gWiz™ High-Expression HBsAg vector (gWiz™ HBsAg)	10 mg	\$850		

Subtotal

Shipping & Handling: Aldevron routinely ships overnight with FedEx, UPS, and DHL. Special shipping needs may require additional charges.

Continental
United States
\$55

CAN/AK/
HI/PR
\$85

International
TBD
approx: \$125

S & H

Purchase Order or Credit Card Number: _____

Expiration Date: _____ Visa MasterCard American Express Discover

I would like to be contacted for payment information. (recommended for email orders)

Total USD

Shipping Info.

I am a new client

Billing Info.

Contact Person	Bill To (name from credit card if used)
Company	Address - line 1
Address - line 1	Address - line 2
Address - line 2	City State Zip/Postal
City State Zip/Postal	Country VAT Exemption # if applicable
Country Telephone Fax	Telephone Fax
email	email